

ADD-ON LABS

Fewer draws, faster results!



If you ordered a BMP, CBC or INR, you can place add-on orders for the tests listed below each one.

BMP (yellow top tube/SST)

Acetaminophen level
 AFP
 Amylase
 ANA
 Bilirubin (fractionated)
 CRP
 C3/C4
 Cardiolipin Ab
 Ceruloplasmin
 Cortisol, serum
 Double-stranded DNA (Anti-DS DNA)
 EBV Ab Panel
 Ethanol
 ENA-10
 Ferritin
 Folate
 Hepatitis Panel
 Hep A Ab Panel
 Hep B Surface Ag/Ab
 Hep B Core Ab Panel
 Hep B e Ag/Ab
 Hep C Ab
 Hepatic function panel (LFTs)
 HIV Ag/Ab
 Iron Panel (Iron, TIBC, transferrin)
 LDH
 Lipase
 Lipid Profile
 Magnesium
 Osmolality, serum
 Parathyroid Hormone (within 24 hrs of draw)
 Phosphorus
 Procalcitonin (within 8 hrs of draw)
 Rheumatoid Factor
 Salicylate level
 T3 (free, total)
 T4 (free, total)
 Thyroglobulin Ab
 TPO Ab
 TSH
 Vitamin B12
 Vitamin D (25-OH)
 Uric acid

CBC (lavender top tube)

BNP (within 8 hrs of draw)
 Carbon Monoxide Level
 CBC (differential)
 ESR
 Hgb A1c
 Retic Count
 Immunosuppressants (Tacrolimus, Cyclosporine)

INR (blue top tube)

aPTT
 Anti-thrombin III Activity
 Anti-thrombin III Ag
 Factor 5
 Factor 8 (along with other factors)
 D-Dimer
 Fibrinogen
 Protein S
 Protein C
 PT/INR
 Add ons must be requested within 4 hrs of collection (PT can be requested up to 24 hr)

QUICK FACTS:

- Chemistry keeps specimens for 4 days (day of collection + 3 days).
- Hematology keeps specimens for ~24hrs only.
- Volume of most tubes is ~3cc and most add on labs only use 0.25 - 0.5cc

Scan for ordering instructions

